


Homogenizer with controller

HG-15A / HG-15D



1 Features

- High speed, high quality and low noise at max. 27,000 rpm
- Speed & power control optimized for the sample
- The dispensing tool made of stainless steel and PTFE
- High efficiency and economical price

2 Specifications

Model	HG-15A	HG-15D
Work range	Max. 2.5 liter	
Speed range	2,000 - 27,000 rpm	
Speed resolution	Stepless	10 rpm
Speed control	Analog phase control	Digital remote control
Usable tool	HT1018, HT1025	
Drive system	Motor / 160W	
Safety device	Over load protector	
Power source	220VAC, 50/60Hz	
Overall (WxBxH)	60x70x230mm	
Weight (Kg)	3	

3 Options

Homogenizer Tool, HT1018, Rotor ϕ 12.7, Stator ϕ 18, 10~1000ml

Homogenizer Tool, HT1025, Rotor ϕ 17, Stator ϕ 25, 50~2500ml

CL100, Clamp for Stirrer or Homogenizer, Cast Aluminum, Grip Capacity ϕ 23mm for Stand/Rod

ST110, Standard Plate Stand for Stirrer and Homogenizer, ϕ 16 x h800 mm Rod, 20 x 35 cm Plate

ST120, Wide Plate Stand for Stirrer and Homogenizer, ϕ 16 x h1000 mm Rod, 30 x 40 cm Plate



HT-1018



HT-1025



CL-100



ST-110

ST-120

Homo Mixer & Homo Disper



HM-1200D

HD-1200D

Rotors



Stators



Dispersing tool

1 Features

- High speed, high quality and low noise at max.
- Speed & power control optimized for the sample
- The tool made of stainless steel and PTFE
- High efficiency and economical price

2 Specifications

Model	HM-1200D	HD-1200D
Work range	500 - 2000ml	500 - 5000ml
Speed range	1000 - 11000 rpm	1000 - 12000 rpm
Speed control	Feed back & Scale	
Speed display	Digital	
Usable tool	Stator :Dia.20/30mm, Rotor :Dia.15/23mm	
Drive system	Motor / 142W	
Safety device	Over load protector	
Power source	220VAC, 50/60Hz	
Overall (WxBxH)	70x180x170mm	
Weight (Kg)	3.2	2.8

3 Options : 1) Dispersing tool 2) Vessel clamp

4 Ordering ; Package / HM-1200D + Stand + Clamp