


VS / IF Series


snowed ice

VS-625NS

flaked ice

IF100**1 Features**

- These units produce flaked or snow type ice for laboratory
- The unit turns off automatically when the storage bin is full and when there is no water
- Model VS-625NS uses a rolling curve to produce the form of snow ice
- Model IF-100 uses a screw turbine system to provide uniform flakes and minimize liquid residue during production

2 Specifications

Model	VS-625NS	IF-100
Ice type	Snow type	Flake type
Production capacity	50 Kg / day	150 Kg / day
Chamber (WxBxH)	300x410x200mm	450x450x500mm
Storage	25 liter	100 liter
Refrigerator	3/8 Hp	1/2 Hp
Water supply	External / Internal supply	
Material	Stainless&Steel+ABS	Stainless & Steel
Door	Front glass door	Front sliding
Safety device	Circuit breaker, Water level & alarm	
Power source	220VAC, 50/60Hz	
Overall (WxBxH)	360x475x615mm	1000x600x1000mm
Weight (Kg)	130	180

Ice Cube Maker

SCI series
**SCI-035****SCI-050 / 090****SCI-120A****1 Features**

- Approved ETL and sanitation, easy volume adjustment controlled by microprocessor.
- Crystal clear ice production and clean storage.
- High density urethane foam insulation and durable stainless steel exterior

2 Specifications

Model	SCI-035	SCI-050	SCI-090	SCI-120A
Capacity / 24 hours (Kg)	32	48	89	125
Storage (Kg)	10	13	45	45
Refrigerator (Hp)	0.37	0.55	0.75	1.0
Power consumption (W)	450	800	1100	1400
Power source	220VAC, 50/60Hz, single phase			
Overall (WxDxH) mm	400x445x770	400x640x950	700x675x1160	750x645x1170
Weight (Kg)	50	58	90	102