



KSU-200



KSU-300 / 400



KSU-600

1 Features

- Study, non-corrosive units handle all models of cleaning tasks
- Operating at 40KHz, these units provide aggressive cleaning powder
- Stainless steel bath have no corner welds and other dirt-cating surfaces
- Models have analog timer, heater and drain valve (KSU-300/400)
- The heat-up temperature is Max.60 °C
- Set the cleaning cycle up to 60 min.



KSU-100

2 Specifications

Model	KSU-100	KSU-200	KSU-300	KSU-400	KSU-600
Capacity	1.8 liter	3.3 liter	6.6 liter	10 liter	20 liter
Bath (WxBxH)	140x135x100mm	240x140x100mm	295x150x150mm	295x240x150mm	498x298x150mm
Temp. range	Amb. to +60 °C				
Timer	0 - 60min				
Wave system	BLT40KHz,60W	BLT40KHz,100W	BLT40KHz,150W	BLT40KHz,200W	BLT40KHz,600W
Drain valve	No		Yes		
Control&display	Control knob & analog				
Material	ABS+Stainless steel		Stainless steel		
Power source	110/220VAC,50/60Hz				
Overall (WxDxH)	160x145x200mm	260x160x200mm	320x160x290mm	320x260x300mm	520x320x320mm
Weight (Kg)	5.5	8	12	18	25

3 Options : Wire Basket

- 1) BTW 300 / 250x100x100H(mm) for KSU-300
- 2) BTW 400 / 250x200x100H(mm) for KSU-400
- 3) BTW 600 / 450x250x100H(mm) for KSU-600

