

Powersonic 212



1 Features

- Automatic water supply and drain by Siphon principle
- Strong cleaning power by industrial BLT transducer
- No need to change water when cleaning
- Short cleaning time Convenient operation plate
- Available to clean multiple pipettes with safety by a basket

2 Specifications

Model	Power Sonic 212
Frequency	28 KHz
Timer	Digital, 0 - 99 min
Vessel	Dia.160 x 685H(mm)
Basket	Dia. 140 x 505H(mm)
Water Flow Rate	12 liter
Cabinet Material	Stainless steel
Vessel Material	Tranparensy acrylic resin
Output	250W
Power source	220VAC,50/60Hz
Overall	320Wx320bx850H(mm)
Weight	15 kg

Auto Spray Hand Sterilizer

CK-003HS



1 Features

- Supply antiseptic solution from well air-tight tank
- The fast natural evaporation become just after spray, and dry a few sec.
- Second generation of contagion won't occurred
- Fume free and an economical price
- Simple supply and drain of antiseptic solution

2 Specifications

Model	CK-003HS
Spraying	1.3 cc/sec
Reservoir	2,000 cc
Material	Stainless steel case & Acryl chamber
Antiseptic solution	Cleancohol (Ethanol 75%)
Control	Power On-off switch, Automatic on-off
Power source	220 VAC,50/60Hzmsingle phase,64W
Overall (WxDxH)	315x300x420mm
Weight (Kg)	8

3 Options : Cleancohol (20 L)