

AND FX series / Japan 

FX-2000i / 3000i



FX-200i / 300i

1 Features

- ▶ Compact B5 size
- ▶ Fast stabilization speed – Just 1 second
- ▶ Multiple weighing units
- ▶ Standard RS-232C interface
- ▶ USB interface – driver installation not necessary
- ▶ Ethernet interface with WinCT-Plus software
- ▶ GLP, GMP, GCP, ISO compliant
- ▶ Statistical calculation function
- ▶ Comparator function
- ▶ Hold function for Animal weighing application
- ▶ Large bright Vacuum Fluorescent Display
- ▶ Built-in Rechargeable Battery
- ▶ Under hook for suspended weighing
- ▶ Standard Breeze break for 1mg resolution models
- ▶ Data memory function
- ▶ Counting function with ACAI
- ▶ Auto Power ON function
- ▶ Auto Power OFF function
- ▶ Digital Calibration Function

2 Specifications

Model	FX-200i	FX-300i	FX-2000i	FX-3000i
Capacity	220 g	320 g	2200 g	3200 g
Gram(g)	0.001 g		0.01 g	
Weighting	g, %, ct, counting			
Linearity	±0.002 g		±0.02 g	
Repeatability	0.001 g		0.01 g	
Stable Time	1 sec			
Display Refresh	5 time / sec, 10 time / sec			
Pan size	Dia. 130 mm		Dia. 150 mm	
Sense Drift	±2 ppm/°C (10°C - 30°C)			
Operating Temp	5°C - 40°C / 41°F - 104°F, Less than 85% RH			
Dimension	190Wx262Bx234H(mm)		190Wx262Bx84H(mm)	
Weight	Approx. 3.0 Kg		Approx. 2.5 Kg	
Power source	220VAC,50/60Hz			
Std. Accessories	Manual, AC adaptor, Fuse, Dust cover			

- 3 Options** : 1) USB interface
 2) Ethernet interface
 3) Built-in Rechargeable Battery
 4) Small Breeze Break
 5) Large Breeze Break