


Colony Counter

HUMAN Lab® DCC-560



1 Features

- This counter features just touch the petri dish with any marker to achieve the counting mechanism
- Counter confirms each mark with an audible signal

2 Specifications

Model	DCC-560
Counting	Integral pressure sensor
Read-out	4 digits LED display
Dish Holder	Petri dish up to 15cm
Illumination	Ring shaped lamp
Count Control	Push button, Needle point probe, Marking pen
Electric power	220VAC, 50/60Hz
Overall (WxDxH)	275x335x175 mm
Weight	5 Kg

3 Accessories ;

- 1) Magnifier (2.5x) & Stand x 1
- 2) Needle point probe x 1
- 3) Marking pen x 1
- 4) Petri dish holder(10cm) x 1

Differential Luokocyte Counter

312DF



1 Features

- 12 Keys programed for counting / identifying items
- 4 function key (% , conti, +/-, clear)
- Range up to 999 for each group and digital display
- Audible sound at 100 counts

2 Specifications

Model	312DF
Numeric key	12 key
Function key	4 of %, conti, ±, clear
Display	Digital, 3 digits
Electric power	110/220VAC, 50/60Hz
Overall (WxDxH)	180x250x85 mm
Weight	1 Kg

Cell Counter

HIP-MC09N



1 Features

- This unit is ideal for cell counting in the lab.
- 12 Keys for counting and function
- Audible sound at 100 counts

2 Specifications

Model	HIP-MC09N
Counting key	9 key
Function key	3 of clear, % & +/-
Display	LCD display
Alarm range	Audible 100 count
Electric power	110-220VAC, 50/60Hz
Overall (WxDxH)	165x235x60mm
Weight (Kg)	0.7