

Desiccator

SK Series

Acryl Auto Desiccator



SK-C015

Acryl Std. Desiccator



SK-C023

Acryl Temp. Desiccator



SK-C016

Acryl Humi. Desiccator



SK-C036

Aluminum Desiccator



SK-C013



SK-C012

Aluminum Auto Desiccator



SK-C003

Vacuum Desiccator



SK-V002



SK-V006

1 Features

- Automatically operate 24 hr-timer and to range 20-50% humidity(SK-C015)
- Designed for const. temp. +5 to 60°C (SK-C016)
- Typical auto.desiccator with hygrometer(SK-C023)
- Designed for humidity chamber (SK-C036)

2 Specifications

| Model | Description |
|---------|-------------------------------------|
| SK-C015 | Auto Desiccator |
| | 600W x 355D x 320H(mm) |
| | Humidity : 20 - 50%, Shelf x 2 |
| SK-C016 | Temperature Desiccator |
| | 620W x 355D x 345H(mm) Shelf x 2 |
| SK-C023 | STD Desiccator |
| | 475W x 365D x 320H(mm) |
| | Shelf x 2, Tray x 1 |
| SK-C036 | Humidity Desiccator |
| | 475W x 365D x 320H(mm) |
| | SUS Net x 5 |

1 Features

- Constructed aluminum frame and PMMA
- Convenience for use shelves, which are adjustable at interval of 50mm height
- Automatically humidity control (SK-C003)
- Designed to allow attachment of variety accessories
- Bigger size and the door is tightly sealed by catch

2 Specifications

| Model | Description |
|---------|-------------------------|
| SK-C013 | 575W x 520D x 770H(mm) |
| | Shelf x 2, Gel tray x 1 |
| SK-C012 | 575W x 520D x 1020H(mm) |
| | Shelf x 4, Gel tray x 1 |
| SK-C003 | 575W x 520D x 1020H(mm) |
| | Shelf x 2, gel tray x1 |
| | 20% ~ Room Humidity |

1 Features

- Useful for testing in the various condition and usable for desiccator alone.
- Designed to reduce pressure up to 1 Torr and utilized for the vacuum storage.

2 Specifications

| Model | Description |
|-------------|------------------------------|
| SK-V002 - A | 260W x 280D x 410H(mm) |
| | 300W x 280D x 510H(mm) |
| | with Shelf x 3, Gel tray x 1 |
| SK-V006 - A | 260W x 180D x 260H(mm) |
| | 260W x 200D x 360H(mm) |
| | 260W x 180D x 460H(mm) |
| | with Shelf x 2, Gel tray x 1 |