



Hahnvapor HS Sereis



HS-2005V-N



HS-2005S-N



HS-5SP

1 Features

- Visual display of rotating speed, temperature of the water bath and vapour
- Excellent in durability of the body and every area of contacting vapour
- PTFE sealing method is adopted for glass contacted area
- Easy for washing and removal of viscous material with evaporating flask of dia. 100mm

2 Specifications

| Model | HS-2005V-N / HS-2005S-N | HS-5SP |
|----------------------|---|----------------------------|
| Rotational Range | 20 ~ 280 rpm | 10 - 180 rpm (Max.280 rpm) |
| Drive System | Brushless DC motor | |
| Lift System | Motorized automatic | |
| Glass Material | Pyrex | |
| Evaporating Capacity | 50 - 3000ml | 5000ml |
| Receiving Capacity | 50 - 3000ml | 3000ml |
| Pressure & Vacuum | Max. 2 bar & 710 mmHg | |
| Bath Capacity | 4.3 liter | |
| Bath Temperature | 20 ~ 180 °C | |
| Control & Display | Digital LED display | |
| Safety Device | Over temp. limiter, Over current protection | |
| Power source | 220VAC, 50/60Hz | |
| Overall (WxDxH)mm | 615x410x880 | 500x700x1600 |
| Weight (Kg) | 24 | 48 |

3 Accessories

- 1) Water Bath / Model : HS-3001
- 2) 1 liter of glass set

4 Options : 1) Electric Aspirator / Model :HS-3000

- 2) Water Circulator / Model : HS-3005N
- 3) Diaphragm Vacuum Pump / HS-3004
- 4) Vacuum Controller / Model : HS-0245