



HS-3000

## Electric Aspirator

Model	HS-3000
Bath Capacity	13 liter (400x240x350mm)
Pump capacity	18 liter/min x 2
Dirve Motor	150W
Material	Polypropylene
Power source	220VAC,50/60Hz
Overall (WxDxH)	450x250x550(mm)
Weight	6.5 Kg



HS-3001

## Heating Water Bath

Model	HS-3001
Capacity	4.3 liter (250Dx120H)mm
Temp. Range	Amb. up to 180 °C
Heater	1000W
Temp. Control	Digital control & display
Material	Bath : Stainelss steel
Power source	220VAC,50/60Hz
Overall (WxDxH)	250x260x190(mm)
Weight	2.8 kg



HS-3005N

## Water Circulator

Model	HS-3005N
Temp. Range	Down to -21 °C
Refrigerator	343W ( 295Kcal / h )
Circulation Pump	Magnet 40W x Max.Flow :19 L/min
Bath Size	6 liter (220Dx210H)mm
Control & Dispaly	Digital
Power source	220VAC,50/60Hz
Overall (WxDxH)	352x410x520(mm)
Weight	26 Kg



HS-3004

## Diaphragm Vacuum Pump

Model	HS-3004
Delivery	20 liter/min
Ultimate vacuum	8 mbar.abs
Operating	1 bar
Drive motor	120W / 0.7A
Operatiom temp.	Amb.+5 to +40 °C
Outlet size	ID 10mm tube
Power souce	220VAC,50/60Hz
Overall (WxDxH)	315x155x210H(mm)
Weight	10 Kg



HS-0245

## Vacuum Controller

Model	HS-0245
Control Range	101 ~ 0.1 PSI ±3% FS
Control & Dispaly	Digital
Power source	220VAC,50/60Hz
Overall (WxDxH)	160x130x95(mm)
Weight	1.7 Kg