

Soxhelt Water Bath

HUMAN Lab® SWB Series



SWB-6D



SWB-4x2

1 Features

- This soxhelt water bath are designed to gently warm samples, evaporate liquids, and melt solids such as agar.
- Concentric ring can be removed from inner cover to accommodate various vessel shapes and sizes
- Four or Six holes model provide bath temperature up to 100°C

2 Specifications

Model	SWB-6	SWB-3x2	SWB-4x2
Capacity	14 liter - 6x1	20 liter-3x2	26 liter-4x2
Chamber (WxBxH)	890x150x100	520x350x110	680x350x110
Temp. range	Ambient to Max. 100°C		
Timer	Digital / 99hr,59min		
Heating	1.5 Kw	1.6 Kw	1.8 Kw
Material	Interior	Stainless steel	
	Exterior	Epoxy coated steel plate	
Control & Display	Microprocessor PID control & digital display		
Safety Device	Over temp. limiter		
Power source	220VAC,50/60Hz		
Overall (WxDxH)	920x180x190	700x420x235	860x420x235
Weight (Kg)	13	16	20


3 Options : Glassware

Extraction Mantle

HUMAN Lab® EAM Series



EAM-9202-03



EAM-9202-06

1 Features

- These extraction mantle are designed for use in distillation, extraction and kjeldahl works
- Each opening has temp. regulator by bimetallic thermocontroller
- The casing is made of stainless steel with excellent chemical resistance
- Moulded silica insulation is achieved the surface temperature are very low

4 Options :

- 1) Glassware set
- 2) With stirring function (EAMS series)
- 3) Digital display (EAMD series)

2 Specifications

Model	Capacity	Heating	Dimension / Weight
EAM-9201-03	100ml x 3	80W x 3	450x240x80 / 4.6Kg
EAM-9202-03	250ml x 3	150W x 3	520x260x95 / 5.8Kg
EAM-9203-03	500ml x 3	260W x 3	580x280x105 / 6.7Kg
EAM-9204-03	1000ml x 3	350W x 3	670x310x115 / 7.9Kg
EAM-9201-06	100ml x 6	80W x 6	950x240x80 / 9.2Kg
EAM-9202-06	250ml x 6	150W x 6	1040x260x95 / 11.6Kg
EAM-9203-06	500ml x 6	260W x 6	1150x280x105 / 13.5Kg
EAM-9204-06	1000ml x 6	350W x 6	1340x310x150 / 15.8Kg

- Max. temp.range : 450°C
- Max. stirring speed : 100 - 1200 rpm (EAMS series only)