


## Systemic Diagram



**Vacuum Pump**

Model	Specification
TRP-06	100 L/min
TRP-12	200 L/min
TRP-24	400 L/min
TRP-36	600 L/min
Vacuum : $1.5 \times 10^{-3}$ Torr	



**Soda Acid Trap**



**Pump Inlet Filter**



**Exhaust Filter**

## Clear Dry Chamber



**General Use**

**Heating plate**

**Stoppering system**

**Complete Flask**

CG3	CH3	CS2	Specification	
Bulk tray unit x 3	Bulk tray unit x 3	Bulk tray unit x 2	CF40 / 40ml	CF600 / 600ml
	Heating plate x 4	Sealing device	CF80 / 80ml	CF900 / 900ml
	Temp.control unit x 1	Stoppering press	CF120 / 120ml	CF1200 / 1200ml
	Temp. sensor x 3		CF150 / 150ml	CF2000 / 2000ml
			CF300 / 300ml	



**Drum Dry Chamber**

**Manifold Dry Chamber**

**Flask Valve**

**Adaptor**

Specification	Flask	Ample	Specification	Specification
FDC-12 : 12 port x 1/2" valve 230Dia.x275H(mm)	FMC4 / 4ea	AMC24 / 24ea	FV 34 / 3/4"	FA-B34 / Bend 3/4
FDC-18 : 18 port x 3/4" valve 300Dia.x320H(mm)	FMC8 / 8ea			FA-S34 / Straight 3/4
	FMC12 / 12ea			AA-34 / Ample 3/4