

## istek DO-30N, 300L, 350L / Korea



DO-30N



DO-300L



DO-350L

## ① Features

## ▶ DO-30N

- Clear LCD display and automatic temp. compensation
- Compensation for salinity and altitude be manual input
- Push button type automatic O<sub>2</sub> calibration - so simple to use
- Most recent calibration data use for direct measurement
- Internal diagnosis for normal operation check
- High sensitive polarographic type electrode
- Two-way RS232 communication
- Automatic power-off for saving battery-life and low battery alarm

## ▶ DO-300L

- Clear LCD display and automatic temp. compensation
- Compensation for salinity and altitude be manual input
- Push button type automatic O<sub>2</sub> calibration - so simple to use
- Most recent calibration data use for direct measurement
- Internal diagnosis for normal operation check
- High sensitive polarographic type electrode
- Current time setting function
- Two-way RS232 communication

## ▶ DO-350L

- Wide graphic LCD display and automatic temp. compensation
- Compensation for salinity and altitude be manual input
- Push button type automatic O<sub>2</sub> calibration - so simple to use
- Most recent calibration data use for direct measurement
- Easily check the status of final calibration
- Internal diagnosis for normal operation check
- High sensitive polarographic type electrode
- Current time setting function
- Alarm sound for high/low limit during DO measurement
- Two-way RS232 communication, and remote controllable

## ② Specifications

| MODEL               |                  | DO-30N                      | DO-300L                                 | DO-350L     |
|---------------------|------------------|-----------------------------|---|-------------|
| Temperature         | Range / Accuracy | -10 to 60 °C ±1 digit       | -10 to 60 °C ±0.4 °C                    |             |
|                     | Resolution       |                             | 0.1 °C                                  |             |
| DO                  | Range / Accuracy |                             | 0.00 to 19.99 mg/L ±0.5%                |             |
|                     | Resolution       |                             | 0.01 / 0.1                              |             |
| O <sub>2</sub>      | Range / Accuracy |                             | 0.0 to 60.0% ±1 digit                   |             |
|                     | Resolution       |                             | 0.1%                                    |             |
| Air Saturation      | Range / Accuracy |                             | 0.0 to 199.9% ±1 digit                  |             |
|                     | Resolution       |                             | 0.1%                                    |             |
| Salinity Correction |                  |                             | 0 to 70 ppt                             |             |
| Temp Compensation   |                  |                             | Auto                                    |             |
| Calibration         |                  |                             | Auto                                    |             |
| Datalog             |                  |                             | 50 Points                               | 100 Points  |
| Print Capability    |                  |                             | Yes                                     |             |
| Display             |                  |                             | Custom LCD                              | Graphic LCD |
| Inputs              |                  | 8 Pin Din                   | BNC, ATC, Power, RS232C                 |             |
| Outputs             |                  | RS232C (Computer / Printer) | Recorder, RS232C (Computer / Printer)   |             |
| Power               |                  |                             | Rechargeable Battery / DC Power Adaptor |             |
| Overall             | (WxDxH)          | 60 x 196 x 30 (mm)          | 270x190x110 (mm)                        |             |
| Weight              |                  | 600 g                       | 720 g                                   |             |