



Vacuum Pump



■ Dimensions

MODEL	A	B	C	D	E	F	G	H	I	J	K	L	M
TRP-6	468	165	253.5	240	120	146	93	140	50	37	28	258	Ø9
TRP-12	468	165	253.5	240	120	146	93	140	50	37	28	258	Ø9
TRP-24	528	205	289.7	310	140	180	141.5	162	70	52.5	34	295	Ø9
TRP-36	556	205	289.7	310	140	180	141.5	187	75	52.5	34	301	Ø9
TRP-48	669.5	249.8	337	440	155	215	227	200	95	45	30	303	Ø14
TRP-60	742.5	279.9	418	550	185	245	300	230	95	45	45	385	Ø14

■ Technical Data

Model		TRP-6	TRP-12	TRP-24	TRP-36	TRP-48	TRP-60
Displacement Speed ℓ/min	60Hz	100	200	400	600	800	1000
	50Hz	80	170	340	500	670	830
Ultimate Pressure Torr	G / ballast C	5×10^{-4}	5×10^{-4}	5×10^{-4}	5×10^{-4}	5×10^{-4}	5×10^{-4}
	G / ballast I	5×10^{-3}	5×10^{-3}	5×10^{-3}	5×10^{-3}	5×10^{-3}	5×10^{-3}
	G / ballast II	1×10^{-2}	1×10^{-2}	1×10^{-2}	1×10^{-2}	1×10^{-2}	1×10^{-2}
Vapour tolerance	Torr	40	40	25	25	40	40
Oil Filing	ℓ	1.1	1.3	1.4	1.8	3.8	5.5
Inlet port	DN	25KF	25KF	25KF	25KF	40KF	40KF
Outlet port	DN	25KF	25KF	25KF	25KF	40KF	40KF
Moter power	Kw	0.4	0.4	0.75	0.75	1.5	2.2
R.P.M of pump	60Hz	1750	1750	1750	1750	1750	1750
	50Hz	1450	1450	1450	1450	1450	1450
Noise level	dB	50	50	52	52	56	56
Weight	Kg	22.5	24	34	38	70	83