


**400 series**

**603**

## 1 Features

- Built-in microprocessor controller
- Finished functions are ultrasonic power, set time, elapsed time, set temperature, current temperature, heater operation status etc.
- Selectable use with ultrasonic strength controller (3 steps)
- Built-in a digital thermometer and timer to set a condition for cleaning
- Easy to read operation status on digital display
- Built-in drain valve

## 2 Specifications

Model	Powersonic 405	Powersonic 410	Powersonic 420	Powersonic 603
Capacity (Liter)	5.7	10	20	3.3
Bath size (WxDxH)	300x155x150mm	300x240x150mm	500x300x150mm	240x140x100mm
Temp. range	Digital, 0 - 50 °C			Digital 0 - +60 °C
Timer	Digital, 9 - 99 min.			
Transducer	Industrial BLT 40 KHz			
Display	400 series - Digital, 500 series - LCD screen			LED screen
Function	Thermostat			Degas
Material	Engineering plastic exterior & stainless steel bath			ABS + SUS
Power source	220VAC, 50/60Hz			
Overall (WxDxH)	500x285x255mm	510x375x275mm	720x450x275mm	325x225x255mm
Weight (Kg)	7	9	14	5.5

## 3 Options : Tray ( Solid / Perforated)

