




LSC Series



1 Features

- Easy to install and movable, measuring kit can detect saturation of the filter
- Safe maintenance of stored bottles is enough to last for a long time due to total elimination of noxious gases, vapour and smell

2 Specifications

Model	LSC-100	LSC-200
Capacity	1 liter x100 bottles	1 liter x 200 bottles
Chamber (WxBxH)	720x450x1800mm	1520x450x1800mm
Blower	Centrifugal, 100m ³ /hr, <50 dB	
Filter	99.99%, < 2000g	
Illumination	FL lamp x 40W	
Shelves	4 ea	8 ea
Materials	Acrylic resin door, Steel sheet exterior	
Control & Display	Microprocessor control & digital display	
Power source	220VAC, 50/60Hz	
Overall (WxDxH)	800x500x2100 mm	1600x500x2100 mm
Weight (Kg)	140	250

Flammable Safety Cabinet



CFS Series



CFS-500 / Yellow



CFS-1100 / Red

1 Features

- Built-in rubber foot and adjustable the level and height
- Door locking device to keep from the stolen
- All air space system is applied to minimize any hazard to outward
- The exhaust hole can be connectable to the duct system
- Door sill protect not to spill the hazard material to outside
- The each shelf is adjustable the height

2 Specifications

Model	CFS-500	CFS-1100
Volume	500 liter	1100 liter
No. of Shelves	3 (standard)	3 (standard)
Shelf Load	Max. 100 Kg	
Material	Epoxy powder coated steel plate	
Color	Red or Yellow	
Overall (WxDxH)	900x530x1160mm	1100x580x1825mm
Weight (Kg)	120	200