

1. Evaporator Flask
2. Evaporator HAHNVAPOR Duct
3. Flange Screw Coupling
4. Vacuum Seal
5. Evaporator Condenser
6. Vacuum Connector
7. Stop - Cock
8. PTFE Hose
9. Receiver Flask Clamp
10. Receiver Flask.
11. Cooling Finger
12. Condenser Cover



Condenser-A

HS-0510



Condenser-B

HS-0511



Vacuum Seal

HS-0170



Joint Clamp

HS-0172




Receiver Flask Ball Clamp

HS-0171




Adapter

HS-0173




⑦ Evaporating flask (24/40)


- 50ml Model : HS-0605
- 100ml Model : HS-0610
- 250ml Model : HS-0625
- 500ml Model : HS-0650
- 1l Model : HS-0700
- 2l Model : HS-0800
- 3l Model : HS-0900



- ① HAHNVAPOR duct (24/40)
 ∅ =312mm Model : HS-0100
 HAHNVAPOR duct (24/40)
 ∅ =200mm Model : HS-0150



- ② Bulbs for flask (24/40)
 ∅ =100mm Model : HS-0200
 ∅ =250mm Model : HS-0201



- ③ Four bulbs for flask (24/40+14/20)
 Model : HS-0300



- ④ Stop cock (19/38)
 ∅ =120mm Model : HS-0400
 ∅ =175mm Model : HS-0401



- ⑤ Connecting adapters (24/40+14/20)
 Model : HS-0500



- ⑥ Vacuum connector
 Model : HS-0160



- ⑧ Receiving flask (S35/20)
 100ml Model : HS-0930
 250ml Model : HS-0940
 500ml Model : HS-0950
 1l Model : HS-0960
 2l Model : HS-0970
 3l Model : HS-0980